## nano F<sup>1</sup>: compressed air & gas filtration

## NMS 0450 450 scfm medical sterile filter


nano F<sup>1</sup> medical sterile filters are specifically designed and manufactured to exceed the requirements for medical gas pipeline systems. Applications include: dental, medical and veterinary.



general characteristics		
rated flow (scfm) (1)	450	
operating limits		
design operating pressure range (psig)	0 to 232	
design operating temperature range (°F)	35 to 176	
housing construction		
materials	pressure die cast aluminum	
finish	E-Coat <sup>™</sup> surface treatment (internal & external) powder paint (external)	
element construction		
media	borosilicate glass microfiber with Nomex <sup>®</sup> protective wrap	
end caps	303 stainless steel	
replacement element (part no)	E 0450 MS	
components		
condensate drain	MDV 50 as standard	
differential pressure indicator / gauge	NDP 1500 as standard	
air circuit		

inlet & outlet connections (NPT) (2)

dimensions and weights	
A (ins)	5.75
B (ins)	2.05
C (ins)	17.01
D (ins)	3.00
E (ins)	21.85
weight (lbs)	9.2



element performance	
DOP efficiency <sup>(3)</sup>	>99.9999%
particle removal	0.01 micron
recommended operating temperature range	35 to 212°F
design operating temperature range	35 to 248°F
maximum autoclave temperature <sup>(4)</sup>	273°F
pressure drop - clean	1.5 psid
pressure drop - recommended replacement	6.0 psid
flow direction through element	outside to inside
maximum steam sterilization cycles	100
maximum element life	6 months or 1000 hours

# namo

nano-purification solutions charlotte, north carolina united states

nano-purification solutions new bethlehem, pennsylvania united states

nano-purification solutions st. catharines, ontario canada

nano-purification solutions gateshead, tyne and wear united kingdom

nano-purification solutions krefeld, germany

tel:	704.897.2182
fax:	704.897.2183
email:	support@n-psi.com
web:	www.n-psi.com

### www.n-psi.com

Experience. Customer. Service.

# Daeo

## F<sup>1</sup>: compressed air & gas filtration

- (1) at 100 psig. For all other pressures, contact support@n-psi.com
- (2) inlet and outlet connections are NPT threaded to ANSI B2.1
- (3) as specified in HTM 02-01 medical gas pipeline systems
- (4) applies to the element only
- (5) technical specifications subject to change without notice. Direct inquiries to support@n-psi.com or contact 704.897.2182

## features

#### housing features

- third party tested and validated in accordance with HTM 02-01, NFPA 99, ISO 7396-1 and AS 2896
- includes element with stainless steel end caps for compatibility with autoclave sterilization
- designed in accordance with ASME pressure vessel standards
- carries CRN (Canadian Registration Numbers) for approved use in every province of Canada
- "dead stop" bowl feature provides full closure without overtightening
- optimized piston seal ensures zero leakage and ease of bowl removal during element replacement
- no tie rod simplifies element change out and reduces access requirements for bowl removal
- comprehensive range of filter mounting accessories
- 10 year housing warranty

## accessories





#### connecting kits

description	fits	part number
NMS connecting kit	NMS 0450	NCK 0700

## wall mounting kits

description	fits	part number
NMS wall mounting kit	NMS 0450	NMK 0700



technical specifications subject to change without notice. Direct inquiries to support@n-psi.com or contact 704.897.2182